

PRODUCT INFORMATION

Uniprot ID	P23510
Common Name	KY-1005;KY1005;SAR-445229
Synonyms	GP34;CD252;OX40L;TXGP1;CD134L;OX-40L;TNLG2B;TNFSF4;Glycoprotein Gp34;CD134 ligand
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	OX40 Ligand
Description	PE-conjugated Anti-OX40 Ligand(amlitelimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

