

PRODUCT INFORMATION

Uniprot ID	P43489
Common Name	BGB-A445;BGB A445
Synonyms	OX40;ACT35;CD134;IMD16;TXGP1L;TNFRSF4
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	OX40
Description	PE-conjugated Anti-OX40(gimistotug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

