

PRODUCT INFORMATION

Uniprot ID	Q96DA0
Common Name	PBP-1510;PBP1510;PBP 1510
Synonyms	EECP;PAUF;JCLN2;HRPE773;PRO1567
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	PAUF
Description	PE-conjugated Anti-PAUF(ulenistamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

