

PRODUCT INFORMATION

Uniprot ID	Q8NBP7
Common Name	PF-04950615, RN-316, RN316
Synonyms	FH3;FHCL3;HCHOLA3;LDLCQ1;NARC-1;NARC1;PC9
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG2
Reactivity	Human
Target	PCSK9
Description	PE-conjugated Anti-PCSK9(bococizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

