

PRODUCT INFORMATION

Uniprot ID	Q15116
Common Name	PSB103;PSB-103;PSB-205 (cocktail with Tuvonralimab);PSB205 (cocktail with Tuvonralimab);QL1706 (cocktail with Tuvonralimab)
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PD-1
Description	PE-conjugated Anti-PD-1(iparomlimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

