

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	Q15116
<b>Common Name</b>	AEX-1188;AEX1188;AEX 1188;INCMGA-0012;INCMGA 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA 12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz
<b>Synonyms</b>	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	PD-1
<b>Description</b>	PE-conjugated Anti-PD-1(retifanlimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

