

PRODUCT INFORMATION

Uniprot ID	Q15116
Common Name	AEX-1188;AEX1188;AEX 1188;INCMGA-0012;INCMGA 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA 12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PD-1
Description	PE-conjugated Anti-PD-1(retifanlimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

