

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	Q15116
<b>Common Name</b>	CHS-007;CHS007;CHS 007;CHS-7;CHS7;CHS 7;JS 001;JS001;JS-001;JS1;JS-1;JS 1;SO-001;SO 001;SO001;SO1;SO-1;SO 1;TAB-001;TAB001;TAB 001;TAB1;TAB-1;TAB 1;Tripleitriumab;Triprizumab;Tuoyi
<b>Synonyms</b>	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	PD-1
<b>Description</b>	PE-conjugated Anti-PD-1(toripalimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

