

PRODUCT INFORMATION

Uniprot ID	Q15116
Common Name	CHS-007;CHS007;CHS 007;CHS-7;CHS7;CHS 7;JS 001;JS001;JS-001;JS1;JS-1;JS 1;SO-001;SO 001;SO001;SO1;SO-1;SO 1;TAB-001;TAB001;TAB 001;TAB1;TAB-1;TAB 1;Tripleitriumab;Triprizumab;Tuoyi
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PD-1
Description	PE-conjugated Anti-PD-1(toripalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

