

PRODUCT INFORMATION

Uniprot ID	Q9NZQ7
Common Name	STI-1014;STI1014;STI 1014;STI-A1014;STI A1014;ZKAB-001;ZKAB001;ZKAB 001;ZKAB-1;ZKAB1;ZKAB 1
Synonyms	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	PD-L1
Description	PE-conjugated Anti-PD-L1(socazolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

