

PRODUCT INFORMATION

Uniprot ID	Q9NZQ7
Common Name	MPDL3280A, RG7446
Synonyms	PD-L1; CD274; B7-H1; PDCD1L1; PDCD1LG1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	PDL1
Description	PE-conjugated Anti-PDL1 (atezolizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

