

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | Q04609 |
| Common Name | J-591, MLN-591, huJ-591 |
| Synonyms | FGCP;FOLH;GCP2;GCPII;mGCP;NAALAD1;NAALAdase;PSM |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Chimeric |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | PSMA |
| Description | PE-conjugated Anti-PSMA(rosopatamab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

