

PRODUCT INFORMATION

Clone ID	6G7
Target	PVRIG
Synonyms	C7orf15;CD112R
Host Species	Rabbit
Description	PE-conjugated Anti-PVRIG antibody(6G7), IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q6DKI7
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Cell surface receptor for NECTIN2. May act as a coinhibitory receptor that suppresses T-cell receptor-mediated signals. Following interaction with NECTIN2, inhibits T-cell proliferation. Competes with CD226 for NECTIN2-binding.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	PE-conjugated

