

PRODUCT INFORMATION

Uniprot ID	Q92854
Common Name	VX15/2503;VX 15;VX-15;VX 15;VX 15/2503;VX-15/2503
Synonyms	A8;GR3;BB18;CD100;COLL4;SEMAJ;coll-4;C9orf164;M-sema-G
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	SEMA4D
Description	PE-conjugated Anti-SEMA4D(pepinemab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

