

PRODUCT INFORMATION

Clone ID	31A10
Target	SEZ6
Synonyms	BSRPC
Host Species	Rabbit
Description	PE-conjugated Anti-SEZ6 antibody(31A10), IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q53EL9
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The protein encoded by this gene is thought to contain five cysteine-rich motifs that are similar to sushi domains, as well as two domains similar to the amino terminal half of the CUB (for complement C1r/C1s, Uegf, Bmp1) domain. Mutations in this gene have been associated with febrile seizures. [provided by RefSeq, Jul 2016]
Usage	Research use only
Conjugate	PE-conjugated

