

PRODUCT INFORMATION

| | |
|---|--|
| Uniprot ID | P37840 |
| Common Name | NEOD-002;NEOD002;NEOD002;NEOD-2;NEOD2;NEOD2;PRX-002;PRX002;PRX 002;PRX-2;PRX2;PRX2;RG 7935;RG7935;RG-7935;RO7046015;RO-7046015;RO7046015 |
| Synonyms | PD1;NACP;PARK1;PARK4 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | SNCA |
| Description | PE-conjugated Anti-SNCA(prasinezumab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

