

PRODUCT INFORMATION

Uniprot ID	P02786
Common Name	AOC1001;AOC-1001
Synonyms	T9;TR;TFR;p90;CD71;TFR1;TRFR;IMD46
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TFRC
Description	PE-conjugated Anti-TFRC(delpacibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

