

PRODUCT INFORMATION

Uniprot ID	Q495A1
Common Name	AB 154;AB154;GS-0154;GS0154
Synonyms	VSIG9;VSTM3;WUCAM
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TIGIT
Description	PE-conjugated Anti-TIGIT(domvanalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

