

PRODUCT INFORMATION

Uniprot ID	Q8TDQ0
Common Name	TSR-022
Synonyms	HAVCR2;TIM3;TIMD3;FLJ14428;KIM3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	TIM3
Description	PE-conjugated Anti-TIM3 (cobolimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

