

PRODUCT INFORMATION

Uniprot ID	Q8TDQ0
Common Name	INCAGN2390;INCAGN02390;INCAGN-2390;INCAGN-02390
Synonyms	TIM3;CD366;KIM-3;SPTCL;TIMD3;Tim-3;TIMD-3;HAVcr-2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TIM3
Description	PE-conjugated Anti-TIM3(verzistobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

