

PRODUCT INFORMATION

Uniprot ID	O00206
Common Name	EB 05;EB05;EB-05;EB 5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101
Synonyms	TOLL;CD284;TLR-4;ARMD10
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TLR4
Description	PE-conjugated Anti-TLR4(paridiprubart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

