

PRODUCT INFORMATION

Uniprot ID	P20333
Common Name	TBC
Synonyms	p75;TBPII;TNFBR;TNFR2;CD120b;TNFR1B;TNFR80;TNF-R75;p75TNFR;TNF-R-II
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TNFRSF1B
Description	PE-conjugated Anti-TNFRSF1B(emunkitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

