

PRODUCT INFORMATION

Uniprot ID	O95150
Common Name	PF 6480605□PF-6480605□RVT-3101,1D1-1.31
Synonyms	TL1;VEGI
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(L117A,L118A,G120A)
Reactivity	Human
Target	TNFSF15
Description	PE-conjugated Anti-TNFSF15(PF-06480605 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

