

PRODUCT INFORMATION

Uniprot ID	O95150
Common Name	PRA-023;PRA 023;PRA023;PRA-23;PRA23;PRA 23
Synonyms	TL1;TL1A;VEGI;TNLG1B;VEGI192A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TNFSF15
Description	PE-conjugated Anti-TNFSF15(tulisokibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

