

PRODUCT INFORMATION

Uniprot ID	A0A0B4J1U6
Common Name	BCD-180;BCD180;BCD 180
Synonyms	TCRBV1S1A1N1;TCRBV9S1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TRBV9
Description	PE-conjugated Anti-TRBV9(seniprutug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

