

PRODUCT INFORMATION

Uniprot ID	Q9NZC2
Common Name	VGL101;VGL-101
Synonyms	PLOSL2;TREM-2;Trem2a;Trem2b;Trem2c
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TREM2
Description	PE-conjugated Anti-TREM2(iluzanebart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

