

PRODUCT INFORMATION

Clone ID	DM74
Target	Trop2
Synonyms	TACSTD2; GA733-1; M1S1; TROP2
Host Species	Rabbit
Description	PE-conjugated Anti-Trop2 antibody(DM74); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	P09758
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.
Usage	Research use only
Conjugate	PE-conjugated

