

PRODUCT INFORMATION

Uniprot ID	Q969D9
Common Name	CSJ-117;CSJ117;NOV-14;NOV14
Synonyms	TSLP
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TSLP
Description	PE-conjugated Anti-TSLP(ecleralimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

