

PRODUCT INFORMATION

Uniprot ID	Q9NP84
Common Name	PDL 192
Synonyms	TNFRSF12A;FGF-inducible 14;FN14;TweakR;CD266
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	TweakR
Description	PE-conjugated Anti-TweakR (enavatuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

