

PRODUCT INFORMATION

Target	GPRC5D
Synonyms	GPRC5D
Description	PE-conjugated Recombinant human GPRC5D protein with C-terminal mouse Fc tag
Delivery	In Stock
Uniprot ID	Q9NZD1
Expression Host	HEK293
Tag	C-Mouse Fc Tag
Molecular Characterization	GPRC5D(Met2-Glu21) mFc(Pro99-Lys330)
Molecular Weight	
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Liquid, PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 4°C in buffer. Do Not Freeze! Ship on blue ice.
Background	The protein encoded by this gene is a member of the G protein-coupled receptor family, however, the specific function of this gene has not yet been determined.
Usage	Research use only
Conjugate	PE-conjugated






Figure 1. Flow cytometry analysis with 1 µg/mL PE-conjugated GPRC5D-mFc protein (PME100559P) on Anti GPRC5D CAR-Jurkat cells (Blue histogram) or negative Jurkat cells (Red histogram).

DIMABIO CONFIDENTIAL

